

# Tristel Duo™

## for Ultrasound

### Microbiological Efficacy Summary

	Test norm	Test type	Organism	Log reduction required	Log reduction achieved	Laboratory	Clean or dirty	Contact time
<b>Sporicidal</b>	EN 13704	Suspension	<i>Bacillus cereus</i>	> 3 log	> 3 log	MGS	Clean 1	30 sec.
			<i>Clostridium difficile</i>	> 3 log	> 5.19 log	HIRL	Clean 1	
				> 3 log	> 5.19 log	HIRL	Dirty 1	
	Babb JR, Bradely CR & Ayliffe GAJ (J. of Hosp. Inf. 1980 1:63-75)	Suspension	<i>Bacillus subtilis var niger</i>	> 3 log	> 6.72 log	HIRL	Clean 1	
			<i>Bacillus subtilis var niger</i>	> 3 log	> 6.66 log	HIRL	Dirty 1	
	EN 14347	Suspension	<i>Bacillus subtilis</i>	> 4 log	> 4.58 log	IRM	N/A	
			<i>Bacillus cereus</i>	> 4 log	> 4.65 log	IRM	N/A	
			<i>Bacillus subtilis</i>	> 4 log	> 4.18 log	MGS	N/A	
<i>Bacillus cereus</i>			> 4 log	> 4.68 log	MGS	N/A		
<b>Mycobactericidal</b>	EN 14563	Carrier	<i>Mycobacterium avium</i>	> 4 log	> 5.2 log	Biotech Germande	Clean 1	30 sec.
			<i>Mycobacterium terrae</i> (Surrogate for <i>Mycobacterium tuberculosis</i> )	> 4 log	> 4.9 log	Biotech Germande	Clean 1	
	EN 14348	Suspension	<i>Mycobacterium avium</i>	> 4 log	> 6.63 log	MGS	Clean 1	
	DGHM Methodology	Carrier	<i>Mycobacterium avium</i>	> 4 log	> 7.76 log	Hygiene Nord GMBH	Dirty 1	
			<i>Mycobacterium terrae</i> (Surrogate for <i>Mycobacterium tuberculosis</i> )	> 4 log	> 7.23 log	Hygiene Nord GMBH	Dirty 1	
	Griffiths et al. Journal of Hospital Infection (1998)	Suspension	<i>Mycobacterium terrae</i> (Surrogate for <i>Mycobacterium tuberculosis</i> )	> 5 log	> 5.82 log	HIRL	Clean 1	
				> 5 log	> 5.69 log	HIRL	Dirty 1	
<b>Virucidal</b>	EN 14476	Suspension	<i>Poliovirus T1</i>	> 4 log	4.26 log	BluScientific	Clean 2	30 sec.
			<i>Adenovirus T5</i>	> 4 log	4.59 log	BluScientific	Clean 2	
			<i>Murine Norovirus</i>	> 4 log	4.17 log	Eurofins	Clean 1	
			<i>Influenza (H1N1)</i>	> 4 log	4.50 log	Eurofins	Dirty 1	
	DVV and RKI	Suspension	<i>Poliovirus T1</i>	> 4 log	> 4.75 log	Mikrolab GMBH	Clean 3	
				> 4 log	> 4.75 log	Mikrolab GMBH	Dirty 2	
			<i>Adenovirus T5</i>	> 4 log	> 6.25 log	Mikrolab GMBH	Clean 3	
				> 4 log	> 6.19 log	Mikrolab GMBH	Dirty 2	
			<i>Polyoma Virus SV40</i> (HPV Surrogate)	> 4 log	> 5.50 log	Mikrolab GMBH	Clean 3	
				> 4 log	> 5.32 log	Mikrolab GMBH	Dirty 2	
			<i>Vaccinia Virus</i>	> 4 log	> 5.00 log	Mikrolab GMBH	Clean 3	
				> 4 log	> 5.38 log	Mikrolab GMBH	Dirty 2	
<b>Fungicidal/Yeasticidal</b>	EN 14562	Carrier	<i>Candida albicans</i>	> 4 log	> 4.4 log	Biotech Germande	Clean 1	30 sec.
			<i>Aspergillus brasiliensis</i>	> 4 log	> 4.36 log	Hygiene Nord GMBH	Clean 1	
			<i>Candida albicans</i>	> 4 log	> 4.08 log	Hygiene Nord GMBH	Clean 1	
	EN 13624	Suspension	<i>Candida albicans</i>	> 4 log	> 4.45 log	MGS	Clean 1	
				> 4 log	> 4.51 log	Hygiene Nord GMBH	Clean 1	
			<i>Aspergillus brasiliensis</i>	> 4 log	> 4.23 log	Hygiene Nord GMBH	Clean 1	
				> 4 log	> 4.33 log	MGS	Clean 1	
			<i>Fusarium solani</i>	> 4 log	> 4.0 log	Huddersfield	Clean 1	
			<i>Aspergillus flavus</i>	> 4 log	> 4.0 log	Huddersfield	Clean 1	
	DGHM Methodology	Carrier	<i>Candida albicans</i>	> 4 log	> 6.47 log	Hygiene Nord GMBH	Dirty 1	

	Test norm	Test type	Organism	Log reduction required	Log reduction achieved	Laboratory	Clean or dirty	Contact time
<b>Bactericidal</b>	EN 14561	Carrier	<i>Staphylococcus aureus</i>	> 5 log	> 5,8 log	Biotech Germande	Clean 1	30 sec.
			<i>Pseudomonas aeruginosa</i>	> 5 log	> 5,5 log	Biotech Germande	Clean 1	
			<i>Enterococcus hirae</i>	> 5 log	> 5,6 log	Biotech Germande	Clean 1	
			<i>Klebsiella pneumoniae</i> Carbapenemase (KPC) Enterobacteriaceae	> 5 log	> 5,57 log	MG5	Clean 1	
			<i>Vancomycin-resistant</i> <i>Enterococcus faecium</i>	> 5 log	> 5,65 log	MG5	Clean 1	
	EN 13727	Suspension	<i>Staphylococcus aureus</i>	> 5 log	> 6,49 log	HIRL	Clean 1	
				> 5 log	> 6,49 log	HIRL	Dirty 1	
			<i>Pseudomonas aeruginosa</i>	> 5 log	> 6,65 log	HIRL	Clean 1	
				> 5 log	> 6,65 log	HIRL	Dirty 1	
			<i>Enterococcus hirae</i>	> 5 log	> 6,59 log	HIRL	Clean 1	
				> 5 log	> 6,59 log	HIRL	Dirty 1	
			<i>Staphylococcus aureus</i>	> 5 log	> 5,36 log	MG5	Clean 1	
			<i>Pseudomonas aeruginosa</i>	> 5 log	> 5,33 log	MG5	Clean 1	
			<i>Enterococcus hirae</i>	> 5 log	> 5,36 log	MG5	Clean 1	
			<i>Staphylococcus aureus</i>	> 5 log	> 5,34 log	MG5	Clean 1	
				> 5 log	> 5,34 log	MG5	Dirty 1	
			<i>Pseudomonas aeruginosa</i>	> 5 log	> 5,47 log	MG5	Clean 1	
				> 5 log	> 5,47 log	MG5	Dirty 1	
			<i>Enterococcus hirae</i>	> 5 log	> 5,43 log	MG5	Clean 1	
				> 5 log	> 5,43 log	MG5	Dirty 1	
			<i>Staphylococcus aureus</i>	> 5 log	> 5,18 log	Hygiene Nord GMBH	Clean 1	
			<i>Pseudomonas aeruginosa</i>	> 5 log	> 5,54 log	Hygiene Nord GMBH	Clean 1	
			<i>Enterococcus hirae</i>	> 5 log	> 5,27 log	Hygiene Nord GMBH	Clean 1	
			Methicillin-resistant <i>Staphylococcus aureus</i> (MRSA)	> 5 log	> 5,34 log	MG5	Clean 1	
				> 5 log	> 5,34 log	MG5	Dirty 1	
			<i>Neisseria gonorrhoeae</i>	> 5 log	> 5,37 log	Huddersfield	Clean 1	
			<i>Gardnerella vaginalis</i>	> 5 log	> 5,39 log	Huddersfield	Clean 1	
			<i>Streptococcus agalactiae</i>	> 5 log	> 5,43 log	Huddersfield	Clean 1	
	DGHM Methodology	Carrier	<i>Pseudomonas aeruginosa</i>	> 5 log	> 7,23 log	Hygiene Nord GMBH	Dirty 1	
			<i>Staphylococcus aureus</i>	> 5 log	> 7,30 log	Hygiene Nord GMBH	Dirty 1	
			<i>Enterococcus hirae</i>	> 5 log	> 7,52 log	Hygiene Nord GMBH	Dirty 1	

## Bespoke Testing

Test type	Organism	Log reduction required	Log reduction achieved	Laboratory	Clean or dirty	Contact time
Suspension	<i>Acanthamoeba castellanii</i> cysts (one of the causative organisms of <i>Acanthamoeba</i> keratitis)	N/A	> 3,0 log	BluTest	Clean 1	30 sec.

### Clean/Dirty Conditions Key

**Clean 1:** 0.3 g/L bovine albumin

**Clean 2:** 0.6 g/L foetal bovine serum + 0.03% bovine serum albumin

**Clean 3:** Aqua bidest

**Dirty 1:** 3.0 g/L bovine albumin + 3 ml/l sheep erythrocytes

**Dirty 2:** 10% Fetal Calf Serum



**Created by:** Tristel Solutions Limited, Lynx Business Park, Cambs, UK, CB8 7NY  
T +44 (0) 1638 721500 - E mail@tristel.com - W www.tristel.com

For Tristel patent information please visit: <http://www.our-patents.info/tristel>

Copyright © Tristel Solutions  
Mkt-Mic-071-1  
August 2016