




2016 | PRODUCT GUIDE

More than 65 products for best infection prevention practices in health care institutions



10 YEAR
ANNIVERSARY
2006-2016

Thank you!





ON THE COVER

04 / Pre-assembled kits of Hygienic Covers® and Hy21® Supports



03 /
Message from the President
Our #1 Priority

Discover the products identified with these symbols



- 06 /** Hy21® Bedpan Support
- 07 /** Hygienic Covers® — Bedpan and Commode Chair
- 08 /** Hy21® Urinal Support for Men and Women
- 09 /** Hygienic Covers® — Urinal
- 10 /** Hy21® Vomit Support
- 11 /** Hygienic Covers® — Vomit
- 12 /** Hygienic Covers® — Toilet Bowl
- 13 /** Super Absorbent Pads
- 15 /** GO VOM® Vomit Bags
- 16 /** HY-SO® Kit
- 17 /** Hanging Scale, Protective Covers and Accessory Plate
- 18 /** GO GLOV® Disposable Toilet Gloves
- 22 /** Hygie Classic® Commode Chair - Steel
- 24 /** Home Care Products
- 26 /** Products for Travel and Outdoor Activities

NEW



- 13** EXTRA absorbent pad
Absorbs more than 1000 ml
- 14** Adjustable Wall Support System
- 16** Hygienic™ Super Absorbent Powder
- 17** Hygienic™ Garbage Bags Made of Recyclable Material
- 18** Wet Wash Glove with Cleansing Solution
- 20** Hygie Alu-Classic® Commode Chair
100% aluminum
- 27** Sustainable Development
Reduce, Reuse, Recycle

Our First Priority: Prevention!

Founded in 2006, Hygie® is an innovative company who has been developing, manufacturing and marketing specialized products to efficiently reduce the spread of bacteria, consequently, reduce the risk for infection in health care. We are well established in Canada, in more than 700 health care facilities and in home care. We now offer a full line of hygiene products for outdoor activities and travel. At Hygie Canada, we understand the acute needs of our customers. We offer a range of revolutionary products in the field of hygiene, designed to make life easier and simpler for all users.

Our company has distinguished itself from the competition with its patented Hy21® technology. Our products consist of Supports and Hygienic Covers® capable of managing body fluids that are partly responsible for nosocomial infections. Our Quebec based company holds patents in 38 countries and we are currently developing the United States, Europe, Asia and Latin America markets while continuing to offer outstanding service and reliability to the Canadian market.

Since 2014, along with our subsidiary company Hy-Industrie, Hygie Canada has been fully vertically integrated, ensuring a Canadian production as well as superior quality assembly. This has allowed us to maintain a firm control over the quality of the products we offer while simultaneously ensuring the most competitive price on the market. In this process of continuous improvement, we focus our efforts and resources on development, marketing and the global ranking of our products for the benefit of our customers around the world.




Éric Tanguay, President



PROHY21® KITS

PRE-ASSEMBLED CONTAINMENT KITS

Pre-assembled kits of Hygienic Covers® and Hy21® Support to help save time and money and more importantly reduce the spread of pathogens that may lead to cross-contamination or the contraction of nosocomial infections. 

Bedpan

200 Hygienic Covers®
10 Hy21® Bedpan Supports
2 Extenders
(for overweight people)
Code: KT-PRO-HY21-B00



10x



10x



2x

Extender

Urinal

200 Hygienic Covers®
20 Hy21® Urinal Supports
Code: KT-PRO-HY21-UL0



10x



20x

Vomit

200 Hygienic Covers®
20 Hy21® Vomit Supports
Code: KT-PRO-HY21-V01



10x



20x



THE HY21® TECHNOLOGY

Designed in Quebec, the Hy21® technology is based on one principle: prevention. This technology is based on a support system made of plastic and Hygienic Covers® made from recycled materials. This makes it safe and easy for health care providers who are called upon to assist patients presenting a variety of infections. Recommended for single-patient use, these products retain bacteria in a closed space therefore helps prevent and limit their proliferation. The range of Hy21® products offer a wide variety of features to better meet the patients' needs both in a hospital environment or in home care services. Here is a list of our various Hy21® Supports:

- Hy21® Bedpan Support
- Hy21® Urinal Support
- Hy21® Vomit Support
- Hy21® GOPEE Support

ISO 13485 Ensuring patient safety

Obtaining ISO 13485 certification proves our commitment to the quality of our medical material. This allows us to demonstrate that our quality management system has been evaluated and found compliant with the regulatory requirements as well as with the needs of our clients.




Patents registered in 38 countries



Questions about the patents? Contact us! 1.866.588.2221 or info@hygie.com

"1 in 25 patients suffer from a nosocomial infection."



BEDPAN SUPPORT

HY21® TECHNOLOGY



Product Code: BA-BASS-HY21-000



Used with Hygienic Covers®, the Hy21® Bedpan Support with handle replaces both the standard and orthopedic bedpan. The innovative Hy21® technology makes manipulation and use easier for all types of patients who require bedpans.

Recommended for single-patient use, the Hy21® Bedpan Support used with the Hygienic Cover® reduces the time needed to dispose of body fluids, as well as washing and disinfecting soiled equipment (one support per patient, per hospital stay).

The supports are made entirely of durable and recyclable polypropylene plastic. They are stackable to optimize storage space.

SPECIFICATIONS

- Support size:
38.5 x 29.6 x 10.6 cm / 15.5 x 11.75 x 4 in
- Units per case: 20



with Hygienic Cover®

HYGIENIC COVERS®

FOR BEDPAN / COMMUNE CHAIR

Product Code: SA-HYGI-CLBI-000




Hygienic Covers® for bedpan and commode chair are the perfect solution to manage and contain body fluids at point of care, while reducing the risks of cross-contamination.

In less than 30 seconds, the super absorbent pad in the Hygienic Cover® will turn, up to 600 ml of body fluids, into gel, thus reducing clean-up and splashing. After use, simply tie and discard in the designated area.


The Hygienic Covers® are made from recycled material and the super absorbent pad from 100% recycled material.

SPECIFICATIONS

- Hygienic Cover® size: 62.2 x 45.7 cm / 24.5 x 18 in
- Unit box size: 20 x 11.4 x 11.4 cm / 8 x 4.5 x 4.5 in
- Units per case: 400 (20 boxes of 20 Hygienic Covers®)



Bedpan extender for overweight people
Placed on the Hy21® Bedpan to increase its height
BA-BASS-HY21-SFO
- Units per case: 10



Hygienic Cover® for orthopedic bedpans
SA-HYGI-CLSB-000
- Units per case: 400
(20 boxes of 20 Covers)

Part of Hygie's containment technology

Absorption
600 ml
30 sec



URINAL SUPPORT

HY21® TECHNOLOGY



Product Code: HY-SUPP-URI3-000



Unique Design
Unisex - Suitable
for men and women



Used with Hygienic Covers®, the Hy21® Urinal Support with handle allows to efficiently manage and keep urine at point of care. It is ideal for patients who are ill or with a loss of mobility. While easy for men to manipulate, this unique design was created to accommodate the female anatomy with a wide opening to avoid spills or leaks.

Recommended for single-patient use, the Hy21® Urinal Support with the Hygienic Cover® reduces the amount of time needed to dispose of body fluids and to wash and disinfect the soiled equipment (one support per patient, per hospital stay).

The handle makes manipulation easy and can be used to hook the urinal support to the bed, keeping the urinal close at hand.

The supports are made entirely of durable and recyclable polypropylene plastic.

SPECIFICATIONS

- Support size:
13 x 9.5 x 24.4 cm / 5.1 x 3.7 x 9.6 in
- Units per case: 20



with Hygienic Cover®

HYGIENIC COVERS®

FOR URINAL SUPPORT

Product Code: SA-HYGI-URI0-000



Hygienic Covers® for Hy21® Urinal Supports is a quick and efficient solution to dispose of urine, while simultaneously avoiding inconveniences, such as cleaning up, splashing and spills, and reduce the risks of cross-contamination.

In less than 30 seconds, the super absorbent pad in the Hygienic Cover® will turn, up to 600 ml of body fluids, into gel. After using, simply tie and discard in the designated area.

Super absorbent pad made of 100% recycled materials.

Optimal use

Adjust the position of bed to ensure that the patient is fully comfortable. The product is more efficient if the patient is in a half-seated position, rather than completely lying down on their back. Lay the super absorbent pad flat in the Hygienic Cover® or fold in two and place it at the bottom of the Hygienic Cover®.

SPECIFICATIONS

- Hygienic Cover® size: 35 x 15.5 cm / 13.8 x 6.1 in
- Unit box size: 25 x 5 x 15 cm / 9.8 x 2 x 6 in
- Units per case: 400 (20 boxes of 20 Hygienic Covers®)



Part of Hygie's containment technology



VOMIT SUPPORT

HY21® TECHNOLOGY



Product code: HY-SUPP-VOMO-000



Used with Hygienic Covers®, the Hy21® Vomit Support with handle is the quick and efficient solution to dispose of body fluids, while simultaneously avoiding all inconveniences such as clean-up, splashing and spills. It is essential for any person with the urge to vomit, whether in an institution, at home or on the road. Particularly recommended for oncology, dialysis, postoperative recovery rooms and also ideal for gastroenteritis.

Recommended for single-patient use, the Hy21® Vomit Support with the Hygienic Cover® reduces the amount of time needed to dispose of body fluids, as well as washing and disinfecting the soiled equipment (one support per patient, per hospital stay).

The supports are made entirely of durable and recyclable polypropylene plastic. They are stackable to optimize storage space.

SPECIFICATIONS

- Support size:
19 x 11.5 x 13 cm / 7.5 x 4.5 x 5.1 in
- Units per case: 20



with Hygienic Cover®

HYGIENIC COVERS®

FOR VOMIT SUPPORT

Product Code: SA-HYGI-VOMO-000



Hygienic Covers® for Hy21® Vomit Support is the quick and efficient solution for anyone with the urge to vomit who wishes to avoid inconveniences such as cleaning up and splashing, thus reducing the risks of cross-contamination.

In less than 30 seconds, the super absorbent pad in the Hygienic Cover® will turn, up to 600 ml of body fluids, into gel, thus reducing clean-up and splashing. After using, simply tie and discard in the designated area.

Super absorbent pad made of 100% recycled materials.

SPECIFICATIONS

- Hygienic Cover® size: 25 x 25 cm / 9.8 x 9.8 in
- Unit box size: 25 x 15 x 5 cm / 9.8 x 6 x 2 in
- Units per case: 400 (20 boxes of 20 Hygienic Covers®)



Part of Hygie's containment technology



HYGIENIC COVERS®

FOR TOILET BOWL

Product Code: SA-HYGI-PCWC-000

The Toilet Bowl Hygienic Cover® with super absorbent pad transforms a traditional toilet into a dry toilet, thus eliminating the need for flushing. Recommended to dispose of urine and stools contaminated by biocides, or radioactive or chemotherapeutic substances. They are ideal when samples need to be collected for analysis, or simply when no water is available.

In less than 30 seconds, the super absorbent pad in the Hygienic Cover® will turn, up to 600 ml of body fluids, into gel, thus reducing clean-up and splashing. After use, simply tie and discard in the designated area.

The Hygienic Covers® are made from recycled material and the super absorbent pad from 100% recycled material.

SPECIFICATIONS

- Hygienic Cover® size: 71 x 45.7 cm / 28 x 18 in
- Units per case: 360 (20 rolls of 18 Hygienic Covers®)



Gelifies fluids to prevent splashes and spills.
600 ml

**LARGE
FORMAT**



SUPER
ABSORBENT
PAD

Essential during:



Residential water shortages



Water shortage in hotel establishments



Military missions



Natural disasters



**Emergency
Preparedness**

**Prevention
Readiness
Intervention
Recovery**

Contact us to find out about these products.
1.866.588.2221 or
info@hygie.com

SUPER & EXTRA ABSORBENT

ABSORBENT PADS IN TWO SIZES - UP TO 1000 ML

Product Code: SA-TAMP-ABSO-020 (9,5 x 5 po)
Product Code: SA-PADS-8X95-000 (9,5 x 8 po)



SUPER ABSORBENT PAD
9.5 x 5 in
Gelifies up to 600 ml

This super absorbent pad will turn up to 600 ml of body fluids, into gel in less than 30 seconds. It is recommended particularly for nosocomial infections, reducing clean-ups and risks of splashing and spills. Practical and efficient, it makes collecting, transportation and disposal much safer. Made from 100% recycled material. Used inside the Hygienic Covers® for Bedpans, Vomit and Urinals.

Case of 400  

SPECIFICATIONS: SA-TAMP-ABSO-400

- Pad size: 22.8 x 12.7 cm / 9.5 x 5 in

- Units per case: 400 (10 packs of 40 units)

Case of 600  

SPECIFICATIONS: SA-TAMP-ABSO-020

- Pad size: 22.8 x 12.7 cm / 9.5 x 5 in

- Units per case: 600 (20 packs of 30 units)



EXTRA ABSORBENT PADS
9.5 x 8 in
Gelifies up to 1000 ml

Larger and more absorbent. The new EXTRA absorbent pad will turn up to 1000 ml of fluids into gel. This new size is recommended particularly for more abundant quantities of fluid, thus reducing clean-ups as well as the risks of splashing and spills. Practical and efficient, it makes collecting, transportation and disposal much safer. Made from 100% recycled material.

Case of 400  

SPECIFICATIONS: SA-PADS-8X95-000

- Pad size: 22.8 x 20.32 cm / 9.5 x 8 in

- Units per case: 400 (10 packs of 40 units)

WALL SUPPORT SYSTEM

ALWAYS WITHIN ARM'S REACH

**NEW
PRODUCT**

**Product Codes: HY-HOLD-PART-001,
HY-HOLD-PART-002, HY-HOLD-PART-003,
HY-HOLD-PART-004, DI-SACS-0000-000**



The wall support system is adjustable based on the customer's needs, is the perfect solution for hanging boxes of Hygienic Covers®, Hy21® Supports such as Bedpan Supports, Urinal Supports, Vomit Supports, and the GO VOM® display. The Support is made from sturdy, cold roll steel with a durable powder coating finish. It contains no porous material, which simplifies clean-up and disinfection. Fastening included.

- 1 SPECIFICATIONS : HY-HOLD-PART-001
Support System for box of Bedpan Hygienic Covers® or 35-unit GO VOM® display
22.20 x 24.75 X 13.33 cm / 8.75 x 9.75 x 5.25 in
- 2 SPECIFICATIONS : HY-HOLD-PART-002
Support System for box of Bedpan Hygienic Covers® and Hy21® Bedpan Support
36.19 x 24.76 X 13.33 cm / 14.25 x 9.75 x 5.25 in
- 3 SPECIFICATIONS : HY-HOLD-PART-003
Support System for box of Urinal and Vomit Hygienic Covers® as well as Hy21® Urinal and Vomit Support
27.30 x 33.65 X 13.97 cm / 10.75 x 13.25 x 5.50 in
- 4 SPECIFICATIONS : HY-HOLD-PART-004
Complete Support System to hold all Hygie products
67.30 x 33.65 X 13.97 cm / 26.5 x 13.25 x 5.5 in
- 5 SPECIFICATIONS : DI-SACS-0000-000
Wall dispenser for Hygienic Covers®
Strong wall bracket made of plastic.
Easy installation and removal
Can hold up to 20 Hygienic Covers®
Fastening not included.
Units per case: 6



GO VOM[®] VOMIT BAGS

GRADUATED WITH ABSORBENT PAD

Product Code: 1-PRBV1 (50 units display)
 Product Code: SA-HYGI-GOVO-000 (35 units display)







NEW
 SIZE
 35 UNITS




Hygienic and practical, the GO VOM[®] vomit bags with absorbent pads are ideal for emergency interventions (hospitals, medical clinics, ambulances), as well as for travel, either by airplane, car, boat, or other. It is also recommended for isolation in the event of nosocomial infection, in oncology and emergency room, as well as during home care services. Convenient for adults and children.

The super absorbent pads in each GO VOM[®] have been designed to, turn up to 300 ml of body fluids, into gel in a few seconds, thus reducing spills. Graduated lines allow for measurements to be taken easily. The pads are made from 100% recycled material.

- 1** SPECIFICATIONS : 1-PRBV1  
 - Transparent bag: 39.5 x 15 cm / 15.5 x 6 in
 - 50-unit display: 18.7 x 14.5 x 20 cm / 7.5 x 5.5 x 8 in
 - Units per case: 400 bags (8 displays of 50 bags)

- 2** SPECIFICATIONS : SA-HYGI-GOVO-000  
 - Transparent bag: 39.5 x 15 cm / 15.5 x 6 in
 - 35-unit display: 20 x 11.5 x 18.5 cm / 8 x 4.5 x 7.25 in
 - Units per case: 420 bags (12 displays of 35 bags)

EASY TO USE



SUPER ABSORBENT POWDER

GELIFIES FLUIDS IN A FEW SECONDS

Product Code: SA-SAPB-450G-012



- Turns liquid into gel in a few seconds
- Capable of absorbing up to 450 times its weight in fluids and soft matter
- Limits the spread of bacteria
- Can be used for multiple purposes such as home, medical or industrial use
- To be applied directly on the fluids

SPECIFICATIONS


- Units per case: 12 bottles of 450g

hygienic[™]

HY-SO[®] Kit

ISOLATION OF BACTERIA

Product Code: HY-HYSO-000-006



The Hy-SO[®] kit is a prevention tool to proactively isolate a patient as soon as he or she presents the first symptoms of a nosocomial infection. The kit prevents the spread of pathogens, which may cause contamination.

The use of Hygienic Covers[®] allows health care personnel to quickly and efficiently finish their clean-up and disinfection tasks.

Hygie's containment technology contributes to the reduction of risks and costs associated with HAI's (healthcare associated infections).

The handle makes it easier for health care personnel to manipulate the support when caring for patients. The edges of the bedpan offer greater stability and improved comfort for patients.

SPECIFICATIONS

- The Hy-SO[®] kit includes: 1 Hy21[®] yellow bedpan support + 6 Hygienic Covers[®]
- Units per case: 6

Product: Scale
BA-KCY0-1000-000

Product: Protective Covers
EN-PRBA-KCY0-100

KCY01 HANGING SCALE

A method to measure urine using Hygienic Covers®

- Maximum capacity: 1.1 kg / 2.7 lbs.
- Scale: 10 g / 0.5 oz.
- 30 cm (11.75 in) in height
- Dual scale graduated (kg and lbs.)
- Units per case: 1



PROTECTIVE COVERS

Transparent protective covers for the Hanging scale.

- No latex
- Units: 100

Product: ♻️

CH-PLAQ-DIST-000 (For standard Hygie Classic® Commode Chair)

CH-PLAQ-DIST-XL0 (For XL Hygie Classic® Commode Chair)

CH-PLAQ-DIST-ALU (For Alu-Classic® Commode Chair)



ACCESSORY PLATE FOR COMMODE CHAIR

- Plastic plate for easy mounting and dismounting.
- To hook on back of a commode chair.
- Can hold up to 20 Hygienic Covers®
- Durable lid with closing latch
- Top opening to allow easy distribution
- Metal support included for disinfecting wipes
- Units per case: 1

Product: ♻️
SA-HYGI-BAG1-400

**NEW
PRODUCT**



HYGIENIC™ GARBAGE BAGS MADE OF RECYCLED MATERIAL

- White garbage bags with closure cords, ultra resistant to any type of waste.
- Made of recycled material
- No super absorbent pad inside
- Can be used for multiple purposes such as home, medical or industrial use
- Thicker plastic film makes garbage bags much more resistant to tear or heavy weight.
- Bag size: 62.2 x 45.7 cm x 24.5 x 18 in
- Units per case: 400 (10 rolls of 40 bags)

WET WASH GLOVES

WITH CLEANSING SOLUTION — SINGLE USE

Product Code: 1-GLOV15



hygienic™

Fragrance free
Latex free

The Hygienic™ Wet Wash Gloves made of polyester are both soft and ultra resistant. Containing no soap or alcohol, they are perfect for bed-ridden patients, incontinent patients, those requiring home care and can even be used during leisure activities. The gloves are soaked in an aloe-based cleansing solution, and are hypoallergenic and dermatologically tested. The cleaning gloves can be warmed up in a microwave oven for greater comfort.

In health care institutions, the bowl-less toilet helps lessen the risks of transmitting pathogens from patients to health care personnel or even to other patients.

SPECIFICATIONS

- Glove size: 23 x 17 cm / 9 x 6.75 in
- Units per case: 210 (15 packs of 14 gloves)

Available 2016

GO GLOV® DISPOSABLE GLOVE

REPLACES WASHCLOTHS — SINGLE USE

Product Code: GA-SEC-01




Ready to use, the GO GLOV® can replace the washcloth for patients' personal hygiene in medical institutions, retirement homes, during home care services and in maternity wards. This single-use product offers a superior level of hygiene and decreases the spread of pathogens, which lead to contamination and HAI's (healthcare associated infections).

Made of very soft, lint-free cloth, it is resistant and hypoallergenic. Its fibres meet the OEKO-TEX® Standard 100 Class 1 requirements, which makes the product safe to use on babies' skin.

SPECIFICATIONS

- GO GLOV® size: 16 x 22 cm / 6.25 x 8.5 in
- Units per case: 1000 (20 packs of 50 units)

A blurred person wearing a blue protective suit and blue shoes is walking away from the camera down a long, brightly lit hospital hallway. The hallway has white walls, a tiled floor, and several doors on the right side. The text is overlaid on the image in a large, bold, italicized font.

*The use of
Hygienic Covers[®]
is the only way
to handle body fluids
that completely respects
all basic, documented
principles.*

*Report from the 2015
Hygienic Cover[®] Pilot Project
Health Care and
Social Services Centre*

HYGIE ALU-CLASSIC®

100% ALUMINUM COMMUNE CHAIR



New technology to prevent infections

At the forefront of the fight against nosocomial infections, Hygie® is proud to present its new aluminum Alu-Classic® Commune Chair, which is devoid of any porous surfaces. Made in Canada, this model was specifically designed for hospitals and long-term care facilities.

- ✓ Aluminum Design
- ✓ Ultra-light, sturdy and stable
- ✓ Adjustable to 3 different heights
- ✓ Seat-floor height of 15 to 20 in
- ✓ Pivoting armrests
- ✓ Non-porous surfaces
- ✓ Rotates on itself
- ✓ Rustproof
- ✓ Easy to clean



	15, 16 and 17-inch Models		18, 19 and 20-inch Models	
Total width	57 cm	22.5 in	57 cm	22.5 in
Total length	76 cm	30.0 in	76 cm	30.0 in
Adjustable seat height	38 cm	15.0 in	47 cm	18.5 in
	40 cm	16.0 in	49.5 cm	19.5 in
	43 cm	17.0 in	52 cm	20.5 in
Lowest height of armrests	92.7 cm	36.5 in	101.6 cm	40 in
Seat size				
Width	48.3 cm	19.0 in	48.3 cm	19.0 in
Depth	44.5 cm	17.5 in	44.5 cm	17.5 in
Maximum load	400 lbs	181.4 kg	400 lbs	181.4 kg



Contact us or visit our website to learn more about the different types of Hygie Alu-Classic® Commune Chairs and their options.



NEW CATALOGUE
Visit hygie.com or contact your Hygie representative

HYGIE ALU-CLASSIC®

100% ALUMINUM COMMUNE CHAIR

Pivoting armrests

The sturdy, secure armrests can be lowered by pressing a button. They also pivot along the seat axis, which frees that area and facilitates patient transfers.

Easy to clean

This commode chair was specially designed to minimize gaps and hard-to-reach spaces, ensuring quick and easy disinfection.

Wheelbase

The wheelbase of the Hygie Alu-Classic® Commode Chair's frame makes it highly stable and secure for the patient, reducing the risk of overturning.

Single-patient seats

The seats reduce the amount of time necessary for the disposal of body fluids and solids. Easily manipulated from the top of the Commode Chair, it simplifies the work of the health care personnel. The use of Hygienic Cover® recommended. (Used with the Alu-Classic® Commode Chair and the Classic® Commode Chair on page 22).

"This commode chair effectively meets the needs of controlling infections in health care institutions because it does not contain any porous material, which facilitates cleaning and disinfection."

France - Nurse specializing in the prevention of infections



Adjustable height

The leg extender allows for optimal, patient-tailored adjustments. Five heights available. This adjustment provides additional comfort for patients and a more comfortable posture for care staff.

Rubber wheels

The 4-inch rubberized wheels offer 360° pivoting action and are equipped with individual manual brakes to fully immobilize the commode chair. The rubber also improves floor grip without leaving scuff marks.

Made in Canada

The lightness of the commode chair's frame is due to its aluminum construction. The slight tilt in the seatback provides better comfort to the patient.



AC-SIEG-HY21-000
Bedpan seat with holes

Bedpan seat with holes and
Hygienic Cover®



1 CH-SIRO-OU SC-001
Round open seat

2 CH-CUVE-HYAL-000
Pail with holes



AC-SIEG-HY21-SFO
Bottomless bedpan seat



AC-SIEG-CUVE-000
Bedpan seat with bottom

HYGIE CLASSIC®

STEEL COMMUNE CHAIR WITH *E-COAT* COATING



Sturdy and Classic

The Hygie Classic® Commode Chair was designed to meet the need of controlling infections in health care institutions. Hygienic Covers® may be used with the Hygie Classic® Commode Chair to limit the spread of pathogens and to prevent contamination of the environment.

The frame of the Hygie Classic® Commode Chair has a greater wheelbase depth than most commode chairs available on the market, which renders it safer and more stable.

Anti-splash PVC bucket (sold separately)

- ✓ Made of steel with *E-COAT* coating
- ✓ Sturdy and stable
- ✓ Adjustable height
- ✓ Easy to clean
- ✓ Pivoting armrests
- ✓ Non-porous surfaces
- ✓ Available in large sizes (24 inches in width)



Pail supports included (rails)



Contact us or visit our website to learn more about the different types of Hygie Classic® Commode Chairs, such as XL sizes and other options.

	16 and 17-inch Models		18, 19 and 20-inch Models	
Total width	54.5 cm	21.5 in	54.5 cm	21.5 in
Total length	76 cm	30.0 in	76 cm	30.0 in
Adjustable seat height	40.5 cm	16.0 in	47 cm	18 in
	43 cm	17.0 in	48.5 cm	19 in
			51 cm	20 in
Lowest height of the armrests	23 cm	9 in	23 cm	9 in
Seat size				
Width	48.3 cm	19.5 in	48.3 cm	19.5 in
Depth	40.5 cm	16 in	40.5 cm	16 in
Maximum load	400 lbs	181.4 kg	400 lbs	181.4 kg

Seats for the Hygie Classic® Commode Chair

- 1 CH-SIRO-OUSC-001
Round, open seat without cover
Recommended for Hygienic Covers®
- 2 CH-SIAL-OUAC-001
Open, elongated seat with cover
- 3 CH-SIAL-FEAC-001
Closed, elongated seat with cover


See Bedpan Seat options on page 21



MANUAL CLEANING

HIGH RISKS OF CONTAMINATION

The disposal of organic waste in toilets and/or the use of spray wands are unacceptable practices due to the risks of contamination and transmission of nosocomial infections.



The use of spray wand

Fluorescent demonstration of splashing and contamination

HOME CARE PRODUCTS TO MEET THE GROWING NEEDS OF THE POPULATION

For a patient or a family member acting as caregiver, Hygie now offers its range of home health care products to cover all specific needs.

This range of products includes the kits and Hygienic Covers® for Bedpans, Urinals and Vomit Supports, as well as the Wet Wash Gloves.

Some products are necessary for recovery; others are intended to make life easier for users. Feel free to contact us with any questions, so that we may better meet your needs and provide you with the necessary explanations on the use of our products.

The products are available in pharmacies, specialty stores and at our online store at www.hygie.com



Hygienic Cover®
with super absorbent
pad





Bedpan, Vomit and Urinal Kits

Simple, easy and efficient
Each kit contains 1 Support and
6 Hygienic Covers®

Code: 1-KBJ04, 1-KVM06, 1-KUL04



**NEW
PRODUCT**
Available
2016

Wet Wash Gloves

A perfect solution for home care
and for busy people

Code: 1-GLOV15



Hygienic Covers®

Simplify home care
Bedpan/Commode Chair, Vomit, Urinal
Box of 10 Hygienic Covers®

Code: 1-EVB10, 1-EVM12, 1-EUR12

PRODUCTS ADAPTED FOR TRAVEL AND OUTDOOR ACTIVITIES

Whether on a car ride with the kids, a hike, a trek, camping in national parks, or on a long trip into the unknown, it is essential for your comfort to have good quality products with you. These new, easy-to-use products will help you avoid unpleasant moments and make your life a lot easier!



For the whole family!

NEW
GO VOM® Vomit bag

Ideal for motion sickness.
Contains a super absorbent pad
Code: 1-GOVOM16
(Case of 16 GO VOM)

NEW
GO PEE® Urinal Kit

To pee wherever you are
Quick, easy and efficient
GO PEE® kit code: 1-ENSURS-5
(Case of 5 kits)

GO PEE® Refill Code: 1-RECUR16
(Case of (16 x 3) GO PEE® refill)

Products available in specialty stores and
at our online store at www.hygie.com



SUSTAINABLE DEVELOPMENT

REDUCE. REUSE. RECYCLE.

Sustainable development is defined as creating products using only methods that reduce negative impact to the environment. This saves energy and natural resources and poses no danger to employees, communities or consumers.



SUPER ABSORBENT PAD 100% RECYCLED MATERIAL

Super absorbent pads are made from 100% recycled material. No chloride, no aluminum, no alcohol and no perfume.

HYGIENIC COVERS® MADE OF RECYCLED MATERIAL

Hygienic Covers® are made from a rich mixture of high quality, recycled polyethylene resins. They are very resistant and offer an increased level of hygiene.

HY21® SUPPORTS 100% RECYCLABLE MATERIAL

All Hy21® Supports from Hygie, such as the Bedpan Support, the Vomit Support and the Urinal Support are made entirely of durable and recyclable polypropylene plastic.

PACKAGING 100% ENVIRONMENTALLY FRIENDLY

Packaging and shipping boxes are made from recycled paper and water-based ink. 100% recyclable



NEW

Wet Wash Glove
For daily body baths
Ideal for health care institutions,
home care and leisure activities
Code: 1-GLOV15
See page 18



**NEW
PRODUCT**
Available
2016

© Hygie Canada inc. All rights reserved. Printed in Canada. ZF-GUD-1-PROD-003



9140 Leduc Blvd.,
Suite 230, Brossard
Quebec J4Y 0L2

Phone: 1-450-444-6777
Toll free: 1-866-588-2221
Fax: 1-450-444-6222

info@hygie.com
www.hygie.com

